
uGun 360

OS-2250 Operation Manual

OSIM
Inspiring well-being

Rev D

1) English

2) Simplified Chinese 简体中文

3) Traditional Chinese 繁體中文

4) Specification 规格说明

OS-2250 uGun 360

Thank you for purchasing the OSIM uGun 360. Please read the operation manual thoroughly for the correct methods of usage, especially taking note of the section on safety.

This manual should be kept on hand for easy reference. It is your guide to a safe and efficient operation.

Contents

Safety	2
Names of parts and functions	5
Operation	7
Product care and storage	9
Troubleshooting	10

IMPORTANT SAFETY INSTRUCTIONS

When using the product, please take special notice of the following basic precautions: Read all instructions before using this product.

DANGER – To reduce the risk of electric shock:

Always unplug this product from the electrical outlet immediately after use and before cleaning.

WARNING – To reduce the risk of burns, fire, electric shock, or injury to persons:

1) Use this product only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
2) The product should never be left unattended when plugged in. Unplug from electrical outlet /power source when the product is not in use, and before putting on or taking off parts.
3) Do not operate the product under a blanket or pillow. Excessive heating may occur and cause fire, burns, electric shock, or injury to persons.
4) Do not use the product if the electrical outlet /power source is loose.
5) Do not connect or disconnect the power plug /USB cable from the electrical outlet /power source with wet hands.
6) When disconnecting the power plug /USB cable, remove it by holding the plug / USB cable head, and not by the cord /cable.
7) Do not operate the product if it has a damaged cord or plug /USB cable, a malfunction, a power failure, or it is dropped into the water. Injury may occur when power is restored unexpectedly. Return the product to OSIM Service Centre for examination and repair.
8) If the power plug /USB cable or power cord is damaged, it must be replaced by OSIM Service Centre or its authorized service agents.
9) Do not carry this product by the power cord /USB cable or use the power cord / USB cable as a handle.
10) Keep the product power cord /USB cable away from heated surfaces.
11) Do not place the power cord /USB cable under the product or any other heavy objects.
12) Do not wind the power cord /USB cable around the product.
13) Do not operate the product with the air openings blocked. Keep the air openings free of lint, hair, and the like.
14) Remove power plug /USB cable from the electrical outlet /power source immediately should liquid be spilled onto the product as this may result in electric shock.
15) Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
16) To stop using the product, turn off all controls and remove the power plug /USB cable from the electrical outlet /power source.
17) Do not drop or insert any object or flammable material into any openings of the product. This may result in electric shock, malfunction, or risk of fire.
18) Children should be supervised to ensure that they do not play with the product.

19) This product can be used by children aged from 8 years and above and persons with reduced physical, sensory, or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Keep children away when the product is in operation.
20) Cleaning and user maintenance of the product shall not be made by children.
21) Do not use or operate the product when your body or hands are wet.
22) Do not use this product outdoors.
23) Do not use the product in an excessively humid or dusty environment as it may result in malfunction.
24) Do not use the product in temperatures higher than 40°C.
25) Do not expose the product to heaters, stoves or direct sunlight.
26) Do not cover the product during operation. This may result in a rise in motor temperature and subsequent malfunction or fire.
27) If you should use the product in a cold room, do not increase the room temperature abruptly. It is recommended that you increase the temperature gradually to a normal level.
28) If the product is brought into a warm environment from a cold place where it was stored, wait for an hour before use as there may be water droplets on mechanical parts due to condensation. Using the product under such condition may result in malfunction.
29) Do not stand on or in product. Use only while seated.
30) Do not dismantle, drop, or subject the product to shock or vibration.
31) Do not rest or place heavy objects on the product as this may result in malfunction.
32) Some noise may be heard from the product while in use. This is due to the structure of the product and is normal.
33) Do not use the product in combination with other therapeutic equipment or electric blanket, etc., as this may result in ineffectiveness or injury.
34) Do not fall asleep when using the product.
35) When using the product for the first time, start with a gentle massage program and low intensity for a short period to understand the features and determine your preferred intensity.
36) Do not apply massage rollers directly on joints of elbows and knees, breasts, or abdomen.
37) Please do not use the product if you are not feeling well or under the influence of alcohol.
38) Each auto massage program is designed to last 15-30 minutes for optimal effectiveness. Ensure that the product rests for 10 minutes before it is being used again as continuous usage over a long period may activate the overheat protection device.
39) For best results, enjoy your massage regularly, while keeping in mind that moderation is key to avoid excessive stimulation and ensure a balanced experience.
40) Avoid food intake for at least an hour before massage /usage.


41) Stop using immediately should you experience discomfort and consult your doctor if you experience any abnormal or severe pain.
42) This product is not intended for self-diagnosis or self-treatment of conditions that should be managed by a qualified health care provider.
43) If you are under any medication or fall into any of the following categories (but not limited to), please consult your health care provider first before use:
 - Heart diseases
 - High, or low blood pressure, deep vein thrombosis (DVT)
 - Pregnant
 - Diabetes
 - Individuals fitted with medical implants e.g. pacemaker, automatic implantable cardioverter defibrillator (AICD)
 - Recently underwent surgical operations
 - Those with any bone conditions like fractures, osteoporosis
44) Do not overcharge otherwise it may damage the battery and affect usage.
45) The instructions for Class III appliances shall state that they must only be supplied with a safety extra-low voltage corresponding to the marking on the appliance.
46) The product must only be serviced or repaired by the OSIM Service Centre. Do not attempt to disassemble or repair the product.
47) This product contains lithium battery that is only replaceable by a skilled technician. Do not attempt to assemble or repair the product. Contact OSIM service centre for any repair.
48) The battery must be removed from the product before it is discarded. The product must be powered off when the battery is removed, and the battery should be safely disposed of. To remove the battery:
 1. Get the OSIM service personnel to remove the battery case.
 2. Remove the fixing screws from the battery case.
 3. Open the battery case and remove the battery.
 4. The battery must be recycled by an OSIM Service Center or disposed of in accordance with local waste management regulations.

SAVE THESE INSTRUCTIONS

This product is intended for household use only
Warranty is subjected to cancellation if the product is used for commercial purposes.

Names of parts and functions

Main unit


*Product images are for illustrative purposes only and are not drawn to scale.

1. Massage head
2. Extended handle design
3. LED indicator (massage intensity and charging status)
4. Massage intensity button
5. On/Off button
6. USB-C port

Names of parts and functions

Accessory

1. Massage heads (Ball, U-Shaped, Bullet, Flat)
2. USB-A to USB-C cable

Note

- Please use only the provided USB cable for charging.
- Product comes with 4 massage heads, do not use other massage heads that are not intended for this product.
- Handle the product and accessories with care to avoid accidental dropping.

Operation

Charging the product

- Product should be fully charged before first use.
- Insert the provided USB cable into the USB port of the product.
- Connect the USB cable to electronic devices such as laptops, USB wall power adapters or portable chargers.
- Do not operate the product while charging or in a wet environment.
- Each full charge (approximately 4 hours) last up to 2.5 hours of usage on average (Usage time varies based on the massage intensity and massage head selected).
- The LED indicator will blink 5 times when the battery is running low. The product will automatically shut down when the battery is flat.
- The LED indicator will show current battery level with running LED lights from bottom to top while charging.
- The LED indicator will remain fully lit when the battery is fully charged.

Start/Stop the massage function

- Attach the desired massage head to the product.
- Press and hold [⏻] button for 1.5 seconds to start the massage at lowest massage intensity (level 1). The LED indicator will show current battery level before it emits a solid white light at level 1.
- Toggle the [◔] button to cycle through the 4 intensity levels.
- Press and hold [⏻] for 1.5 seconds at any time to power off the product.

Note

- Do not use the product to massage on one spot for more than 30 seconds.

Operation

Recommended usage area of massage heads

- Ball massage head
 Suitable for arms, waist, back, butt, thighs and calves.
- U-shaped massage head
 Suitable for neck, spine, heel and leg.
- Bullet massage head
 Suitable for impacting deep tissues such as meridians, palms and soles.
- Flat massage head
 Suitable for relaxing and massaging large muscles such as the latissimus dorsi.

⚠ Caution

- Recommend to start with low intensity and slowly increase to max intensity after the muscle "loosen" up.
- If you have or have had any medical condition, please seek advice from doctor or other qualified health professional before using the product.
- Stop operation immediately if experience any discomfort.
- It is recommended to wear proper workout attire when using product.

Product care and storage

Product placement and storage

Environmental factors can strongly affect the product physically and mechanically. Do not expose the product to direct sunlight or in places with extreme heat or air-conditioning vents. These can cause the handle grip to dry up, crack or the color to fade.

If the product is not to be used for a long period of time, cover it with a dust cover.

Clean the product regularly using a dry cloth

Before cleaning, ensure that the main switch is turned off and the cable is disconnected from the product.

The massage head of the product is customarily exposed to the natural body oils of users. The whole unit can be cleaned thoroughly with a piece of soft, dry cloth.

The whole unit can be cleaned thoroughly with a piece of soft, dry cloth. The common micro-fibre cleaning materials or gloves are very well suited for this routine maintenance process. Such materials and gloves are available for sale in most DIY shops.

Do NOT use wet elements for cleaning

The product is made up of electrical components and mechanical parts, which are all very sensitive to moisture. Under no circumstances should they be exposed to any form of water or moisture.

Do NOT use cleaning solvents

Under no circumstances should you use any cleaning agents, like benzene or thinner to clean the product. This may cause the color to fade.

Do NOT attempt to repair the product

Under no circumstances should you attempt to repair the product, for safety reason.

If you suspect a malfunction, always call the OSIM Service Center for further assistance and information.

Troubleshooting

If you encounter problems operating the product, please check the following points for rectification. If the problems persist, please switch off the product, and contact the OSIM Service Center* for servicing.

Do not disassemble or attempt to repair the product yourself. Doing so will void the warranty and OSIM will not be responsible for any subsequent consequences.

Problem	Possible Cause / Action
• The product will not start.	• Ensure that you press and hold [⏻] button for 1.5 seconds. • The battery level is low. Recharge the battery before attempting to use the product again. • The overcurrent protection device may have been activated. Turn off the product and allow it to cool down for at least 30 minutes before using it again.
• The massage stops midway.	• The auto timer of 15 minutes may have lapsed. • The battery level is low. Recharge the battery before attempting to use the product again. • The force applied on muscle is too strong and cause the safety features to activate. Reactivate product to assure function. • The overheat protection device may have been activated. Turn off the main switch and allow the product to cool down for at least 30 minutes before using it again.

Troubleshooting

Problem	Possible Cause / Action
• The product gets hot in hand. • The power cord and plug are abnormally hot.	• Stop using product immediately. The overcurrent protection device may have been activated. • Turn off the product and allow it to cool down for at least 30 minutes before touching it again. • If the same problem persists, contact OSIM Service Center* immediately.

* Please visit us at **www.OSIM.com** for information on the nearest service center. Product warranty is valid only in the country of purchase.

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備考 1. 〝超出 0.1 wt %〞及〝超出 0.01 wt %〞係指限用物質之百分比含量超出百分比含量基準值。
Note 1. "Exceeding 0.1 wt %" and "exceeding 0.01 wt %" indicate that the percentage content of the restricted substance exceeds the reference percentage value of presence condition.

備考 2. 〝○〞係指該項限用物質之百分比含量未超出百分比含量基準值。
Note 2. "○" indicates that the percentage content of the restricted substance does not exceed the percentage of reference value of presence.

備考 3. 〝-〞係指該項限用物質為排除項目。
Note 3. The "－" indicates that the restricted substance corresponds to the exemption.

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Designed & Developed in Singapore by OSIM International Pte. Ltd.
© 2025 All rights reserved. www.OSIM.com
OSIM is Subsidiary of V3 Brands