
uGun Mini

OS-2253 Operation Manual

OSIM
Inspiring well-being

Rev D

1) English

2) Simplified Chinese 简体中文

3) Traditional Chinese 繁體中文

4) Specification 规格说明

OS-2253 uGun Mini

ENGLISH

Thank you for purchasing the OSIM uGun Mini. Please read this operation manual thoroughly for the correct methods of usage, especially taking note of the section on safety.

This manual should be kept on hand for easy reference. It is your guide to a safe and efficient operation.

Contents

Safety	2
Names of parts and functions	7
Operation	9
Product care and storage	13
Troubleshooting	15

IMPORTANT SAFETY INSTRUCTIONS

When using the product, please take special notice of the following basic precautions:
Read all instructions before using this product.

DANGER – To reduce the risk of electric shock:

Always unplug this product from the electrical outlet immediately after use and before cleaning.

WARNING – To reduce the risk of burns, fire, electric shock, or injury to persons:

1) Use this product only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
2) The product should never be left unattended when plugged in. Unplug from electrical outlet /power source when the product is not in use, and before putting on or taking off parts.
3) Do not operate the product under a blanket or pillow. Excessive heating may occur and cause fire, burns, electric shock, or injury to persons.
4) Do not use the product if the electrical outlet /power source is loose.
5) Do not connect or disconnect the power plug /USB cable from the electrical outlet /power source with wet hands.
6) When disconnecting the power plug /USB cable, remove it by holding the plug /USB cable head, and not by the cord / cable.
7) Do not operate the product if it has a damaged cord or plug /USB cable, a malfunction, a power failure, or it is dropped into the water. Injury may occur when power is restored unexpectedly. Return the product to OSIM Service Centre for examination and repair.
8) If the power plug /USB cable or power cord is damaged, it must be replaced by OSIM Service Centre or its authorized service agents.
9) Do not carry this product by the power cord /USB cable or use the power cord /USB cable as a handle.

10) Keep the product power cord /USB cable away from heated surfaces.
11) Do not place the power cord /USB cable under the product or any other heavy objects.
12) Do not wind the power cord /USB cable around the product.
13) Do not operate the product with the air openings blocked. Keep the air openings free of lint, hair, and the like.
14) Remove power plug /USB cable from the electrical outlet / power source immediately should liquid be spilled onto the product as this may result in electric shock.
15) Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
16) To stop using the product, turn off all controls and remove the power plug /USB cable from the electrical outlet /power source.
17) Do not drop or insert any object or flammable material into any openings of the product. This may result in electric shock, malfunction, or risk of fire.
18) Children should be supervised to ensure that they do not play with the product.
19) This product can be used by children aged from 8 years and above and persons with reduced physical, sensory, or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Keep children away when the product is in operation.
20) Cleaning and user maintenance of the product shall not be made by children.
21) Do not use or operate the product when your body or hands are wet.
22) Do not use this product outdoors.
23) Do not use the product in an excessively humid or dusty environment as it may result in malfunction.
24) Do not use the product in a damp, poorly ventilated, or high-temperature environment.
25) Do not expose the product to heaters, stoves or direct sunlight.

26) Do not cover the product during operation. This may result in a rise in motor temperature and subsequent malfunction or fire.
27) If you should use the product in a cold room, do not increase the room temperature abruptly. It is recommended that you increase the temperature gradually to a normal level.
28) If the product is brought into a warm environment from a cold place where it was stored, wait for an hour before use as there may be water droplets on mechanical parts due to condensation. Using the product under such condition may result in malfunction.
29) Do not dismantle, drop, or subject the product to shock or vibration.
30) Do not rest or place heavy objects on the product as this may result in malfunction.
31) Some noise may be heard from the product while in use. This is due to the structure of the product and is normal.
32) Do not use the product in combination with other therapeutic equipment or electric blanket, etc., as this may result in ineffectiveness or injury.
33) Do not fall asleep when using the product.
34) When using the product for the first time, start with a gentle massage program and low intensity for a short period to understand the features and determine your preferred intensity.
35) Please do not use the product if you are not feeling well or under the influence of alcohol.
36) Each auto massage program is designed to last 15 minutes for optimal effectiveness. Ensure that the product rests for 10 minutes before it is being used again as continuous usage over a long period may activate the overheat protection device.
37) For best results, enjoy your massage regularly, while keeping in mind that moderation is key to avoid excessive stimulation and ensure a balanced experience.
38) Avoid food intake for at least an hour before massage / usage.


39) Stop using immediately should you experience discomfort and consult your doctor if you experience any abnormal or severe pain.
40) This product is not intended for self-diagnosis or self-treatment of conditions that should be managed by a qualified health care provider.
41) If you are under any medication or fall into any of the following categories (but not limited to), please consult your health care provider first before use:
 - Heart diseases
 - High, or low blood pressure, deep vein thrombosis (DVT)
 - Pregnant
 - Diabetes
 - Individuals fitted with medical implants e.g. pacemaker, automatic implantable cardioverter defibrillator (AICD)
 - Recently underwent surgical operations
 - Those with any bone conditions like fractures, osteoporosis
42) Do not overcharge otherwise it may damage the battery and affect usage.
43) This appliance shall be supplied only with a safety extra-low voltage corresponding to the rating marked on the appliance.
44) The product must only be serviced or repaired by the OSIM Service Centre. Do not attempt to disassemble or repair the product.
45) This product contains lithium battery that is only replaceable by a skilled technician. Do not attempt to assemble or repair the product. Contact OSIM Service Centre for any repair.

46) The battery must be removed from the product before it is discarded. The product must be powered off when the battery is removed, and the battery should be safely disposed of. To remove the battery:
 1. Get the OSIM service personnel to remove the battery case.
 2. Remove the fixing screws from the battery case.
 3. Open the battery case and remove the battery.
 4. The battery must be recycled by an OSIM Service Centre or disposed of in accordance with local waste management regulations.

SAVE THESE INSTRUCTIONS

This product is intended for household use only
Warranty is subjected to cancellation if the product is used for commercial purposes.


Name of parts and functions

*Product images are for illustrative purposes only and are not drawn to scale.

1. Flat smooth glide head
2. Massage gun
3. Massage intensity/charging LED indicator
4. Power/massage intensity button
5. USB-C port

Name of parts and functions

①

②

1. Massage head silicone attachments
 - U-shape target head
 - Round comfort head
 - Cone precision head
2. USB-A to USB-C cable

Note

- Please use only the provided USB cable for charging.
- Product comes with 3 massage head attachments, do not use other massage heads that are not intended for this product.
- Handle the product and accessories with care to avoid accidental dropping.

Operation

Charging the product

- Product should be fully charged before first use.
- Insert the provided USB cable into the USB port of the product.
- Connect the USB cable to electronic devices such as laptops, USB wall power adapters or portable chargers.
- Do not operate the product while charging or in a wet environment.
- Product will not work when it is charging.
- Each full charge (approximately 2 hours) provides up to 1.5 hours of usage. (Usage time may vary on massage intensity and environmental temperature).

LED Charging Indicators

- When the battery is low, the LED indicator will flash according to its current massage intensity.
- When the battery is fully depleted, the product automatically shuts down.
- While charging, the LED indicator will flash white.
- When fully charged, the LED indicator will turn solid white.


Note

- Charge the product every 5 months even if it is not in use, to ensure the battery functions properly.

Operation

Fitting the Massage Head

- When attaching the massage head attachment, slight tilt it and align the edges before inserting it fully onto the massage head.


⚠ Caution

- Do not attempt to remove or change the massage head while the massage gun is operating. Always switch off the device first.

Operation


Start/Stop the massage function

- Press and hold the power button to turn on the product and it will be at massage intensity level 1.
- Press the power button again to cycle through the different levels of massage intensity.

- Auto timer is 15 minutes.
- Press and hold the power button to turn off the product.

Recommended muscle group area


Note

- Do not use the product to massage on one spot for more than 30 seconds.

Operation

⚠ Caution

- Do not massage on the any muscle group for more than 10 minutes.
- Recommend to start with low intensity and slowly increase to max intensity after the muscle group “loosen” up.
- Do not use product on neck, head or any hard areas of the body. Do not use on bones. Use the product only on soft tissue of body.
- If you have or have had any medical condition, please seek advice from your doctor or other qualified health professional before using the product.
- Stop operation immediately if you experience any discomfort.
- It is recommended that you wear proper workout attire when using product.
- To avoid injury, do not touch the parts as shown below while the product is in operation.

Product care and storage

Product placement and storage

Environmental factors can strongly affect the product physically and mechanically.

Do not expose the product to direct sunlight or in places with extreme heat or air-conditioning vents. These can cause the material to dry up, crack or the color to fade.

If the product is not to be used for a long period of time, cover it with a dust cover.

Clean the product regularly using a dry cloth

Before cleaning, ensure that the main switch is turned off and the cable is disconnected from the product.

The massage head of the product is customarily exposed to the natural body oils of users.

The product can be cleaned regularly with soft, dry cloth to wipe off dust and dirt.

The common micro-fibre cleaning materials or gloves are very well suited for this routine maintenance process. Such materials and gloves are available for sale in most DIY shops.

Do NOT use wet elements for cleaning

The product is made up of electrical components and mechanical parts, which are all very sensitive to moisture. Under no circumstances should they be exposed to any form of water or moisture.

Do NOT use cleaning agents

Under no circumstances should you use any cleaning agents, like benzene or thinner to clean the product.

This may cause the color to fade.

Product care and storage

Do NOT attempt to repair the product

Under no circumstances should you attempt to repair the product, for safety reason.

If you suspect a malfunction, always call the OSIM Service Centre for further assistance and information.

Troubleshooting

If you encounter problems operating the product, please check the following points for rectification. If the problems persist, please switch off the product, unplug the USB cable, and contact the OSIM Service Centre* for servicing.

Do not disassemble or attempt to repair the product yourself. Doing so will void the warranty and OSIM will not be responsible for any subsequent consequences.

Problem	Possible causes / actions
• The product will not start.	• Ensure that you press and hold the power button for 2 seconds. • The battery level is low. Recharge the battery before attempting to use the product again. • Ensure that the product is not in charging mode.
• The massage stops midway.	• The auto timer function of 15 minutes may have lapsed. • The battery level is low. Recharge the battery before attempting to use the product again. • The force applied on muscle is too strong and cause the safety features to activate. Reactivate product to assure function. • If the product stops due to overheating, wait at least 30 minutes before attempting to use it again.

Troubleshooting

Problem	Possible causes / actions
• The product gets hot in hand.	• Stop using product immediately. The overcurrent protection device may have been activated. • Turn off the product and allow it to cool down for at least 30 minutes before touching it again. • If the same problem persists, contact OSIM Service Centre* immediately.

* Please visit us at **www.OSIM.com** for information on the nearest service centre. Product warranty is valid only in the country of purchase.

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備考 1.〝超出 0.1 wt %〞及〝超出 0.01 wt %〞係指限用物質之百分比含量超出百分比含量基準值。
Note 1. "Exceeding 0.1 wt %" and "exceeding 0.01 wt %" indicate that the percentage content of the restricted substance exceeds the reference percentage value of presence condition.

備考 2.〝○〞係指該項限用物質之百分比含量未超出百分比含量基準值。
Note 2. "○" indicates that the percentage content of the restricted substance does not exceed the percentage of reference value of presence.

備考 3.〝 - 〞係指該項限用物質為排除項目。
Note 3. The "－" indicates that the restricted substance corresponds to the exemption.

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Designed & Developed in Singapore by OSIM International Pte. Ltd.
© 2025 All rights reserved. www.OSIM.com
OSIM is Subsidiary of V3 Brands